

Back to Basics-: Supervisors Guidelines

Jan Pollard
Chief Audiologist
Sonic

www.sonici.com

 **SONIC**
Everyday Sounds Better

➔ Knowledge Base

- **Developing an understanding of theory and its application is critical- get the foundations right**
- **It is possible to just “press the buttons” but great outcomes need more- thinking clinicians!**
- **Time under supervision is an opportunity to develop the new clinician- they do not get the opportunity after that**

➔ Supervisors

- Usually haven't had training in how to be a supervisor
- Train the way that they were taught
- Train on what they do...but is that Best Practice?
- Sometimes forget the basic theory because it has been a long time
- Requires continuous revision and review
- **START OFF** as you intend to do long term

➔ Teaching and Training

- Involves you talking to them but then repeating back to you what they heard
- Adult learning needs hands on training as well. They needs to show you that they can do it
- Have a clarity of learning objectives, expectations and competency sign-off

➔ Audiometry- basics

- Masking (of course)- inserts & BC masking
- **BC placement**
- **Inserts vs Headphones- Why/When/How**
- **Cochlear dead regions**

→ **Inserts vs Headphones: Why/When/How**

- **Avoid canal collapse**
- **Larger inter-aural attenuation- less masking**
- **Need to be properly inserted**

- **Need to ensure Fitting software option is correct**
- **Need to ensure Targets set up in REM**



Options

View | Audio | Live Speech | REM | MHA/HLS | Tymp

Live Speech

- Show Fall Sounds
- Show Peak Sounds

Analyzer

Octave Bands: 1/3

- Interpolate Curves

Show Targets:

NAL-NL2 Setup

- Show MCL
- Show Average UCL

Percentiles

Default Stimulus: Noise

Default Level: 65

NAL-NL2 Setup

Client

Client Age: 50 years, 0 months

Gender: Unknown Language: Non-Tonal HI Experience: Experienced

Transducer: Insert Tip

- Insert Tip
- Insert Mold
- Supra-Aural
- Loudspeaker 0°
- Loudspeaker 45°

Hearing Aid

Circuit Type: WDRC

Style: BTE

Use BC

Correct by Client REUG

Binaural

Acoustics

Venting Size: Tight

Tube: None

Fitting Depth: Standard

Compression Speed: Fast

Targets

New Level: Levels: 65dB

Add ->

<- Remove

Speech Targets

Show Targets For: All Levels

Show Targets As: Speech-o-Gram & Line

More...


About OK Cancel

➔ **RED flags for Cochlea Dead Zones**

- * hearing loss $> 70\text{dB}$
- * a steeply sloping, pronounced high-frequency hearing loss
- * a reverse curve hearing loss
- * a cookie bite audiogram
- * a patient report of a non-tonal, scratchy type of sound quality
- * a word-understanding score that is much worse than the audiogram suggests it should be.

➔ Travelling Wave Overlap


'Off frequency' hearing



Do these sound like tones or are they scratchy noise?

➔ Reverse Slope Loss- hard to fit

Do these
sound like
tones or are
they scratchy
noise?



Beware the Upward spread of
Masking



➔ Tympanometry and reflexes

- Otoscopy first
- Tymp- Classification and interpretation
- **Expected effects on audiogram- need to look for consistency in the test battery**
- Reflexes- ipsi? Contra? When and why
- Interpretation


➔ **Speech Discrimination**

- **What levels to test and why- steeply sloping losses, roll-over, correlation to PTA**
- **Maximum scores re Hearing Loss**

➔ The Performance Intensity function



➔ Typical PI functions



NORMAL



➔ MS scores based on 3FAHL

| | | |
|----------------|----------------------|-------------|
| Conductive | | 100% |
| Sensori-neural | Av. Loss up to 40 dB | 90- 100% |
| | 41- 50 | 80% |
| | 51- 60 | 70% |
| | 61- 70 | 60% |
| | 71+ | 50% or less |

* OHS
610/810

➔ **COSI- SMART goals: WHY?**

- **Most important part of the consultation, usually the worst done**
- **SMART goals**
- **What does the CLIENT want to achieve- not problems**


➔ Hearing Aids

- Channels vs Bands
- Loudness vs Clarity
- Compression
- Gain/Frequency Response – prescription, venting, NR, FBC, Adaptive
- Real Ear Verification- Open Fit

➔ Channels vs Bands-always confused

BAND

Frequency region where gain adjustment is made



CHANNEL


Frequency region where the same signal processing takes place

➔ Compression

EXPANSION/SQUELCH- increasing gain with increasing input

LINEAR- a 1:1 compression ratio

WDRC- decreasing gain with increasing input




➔ Linear amplification

- The gain stays the same as the input increases
- **Eg Gain = 30 dB**
 - Input 70 dB Output 100 dB
 - Input 80 dB Output 110 dB
 - Input 90 dB Output 120 dB

➔ Input- Output Curves

- A graph that shows us how much sounds are turned up



$$\text{Output} - \text{Input} = \text{GAIN}$$


➔ **What is Compression**

- the range of input sound intensities is “squashed” into a smaller range of output intensities
- Normal hearing people have a wide Dynamic Range, and as hearing loss gets worse this becomes more and more reduced
- e.g. a range of input intensities from 0 to 120 dB SPL may be compressed into an output range of 50 to 100 dB SPL
- **Wide Dynamic Range Compression WDRC** is compression that occurs for a wide range of inputs

→ Compression- 1 dB impact

Normal

Impaired





Linear



Peak Clipping



Compression



If too compressed you get distortion



➔ Compression Ratio- Gain change

It tells us how squashed up we are making the input

Change in Input level

= Compression ratio

Change in Output level

120 dB


40 dB

$$120 / 40 = 3:1$$

So, every 3dB squashes into 1 dB

➔ Compression Kneepoint

- The input level after which compression starts- this may just be moving from one compression ratio to another
- It is the bend in the Input/Output curve



Below a certain signal intensity the amplifier behaves linearly, above this intensity the compression operates

➔ Definitions

- **Gain**

- A difference measure
eg REUG, REAG

- **Response**

- An output (total) measure eg REUR, REAR



➔ How much Gain?

- Lybarger (1944, cited in Byrne 1983) developed the "half-gain rule", in which recommended gain levels equalled half the amount of hearing loss at each frequency. This procedure resulted in more comfortable listening levels than did the "mirroring" method, and has since been used as the basis for many of the prescriptive fitting formulas that were introduced in the late 1970s and early 1980s.

➔ NAL-NL1/2

- The aim of the original NAL prescriptive formula, which has remained constant through a number of revisions, was to maximize speech intelligibility at a listening level that is preferred by the hearing aid user
- Intelligibility is assumed to be maximized when all bands of speech are perceived by the listener to have the same loudness; that is, when the goal of loudness equalization for speech bands has been achieved

➔ **Things that affect GAIN**

- **Venting, tubing size, dome- need to set this properly in the Fitting software**
- **Directivity**
- **Feedback cancellation FBC**
- **Noise Reduction**
- **Automatic Environmental programs**

SONIC | pep.

Thank You



www.sonici.com

 **SONIC**
Everyday Sounds Better